



# Pediatric Crash Cart Refill Kit

#LC017911

LC017901 - Drawer 1 - Simulated Medications
(2) Adenosine Syringe
(1) Atropine Syringe
(1) Amiodarone Vial
(1) Calcium Chloride Vial
(1) Dextrose 50% Vial
(1) Dextrose 50% Syringe
(2) Epinephrine 1:10000 Syringe
(1) Epi-Pen
(1) Lidocaine 2% Syringe
(1) Naloxone (Narcan) 0.4 mg/ml Vial
(1) Naloxone (Narcan) 1 mg/ml Syringe
(3) Sodium Bicarbonate 4.2% Syringe
(1) Broselow tape

The 9 Drawer Pediatric Cart Refill Kit comes with all of the supplies necessary to restock drawers 1 through 8 in a standard Pediatric/Braslow cart. Each drawer contains supplies designed for a different patient size and weight, making the proper items much easier to find. Explore our predesigned kits or design your own!



\*Cart Not Included in Refill Kit

LC017902 - Drawer 2 - Patient Weight: 0.5 kg - 8 kg	
VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 3.0 mm
(2) IV Catheter: 22 Gauge	(1) Cuffed Endotracheal Tube: 4.0 mm
(2) IV Catheter: 24 Gauge	(1) Small Stylet: 6 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 6 FR
(2) IV Extension Set: Male to Female	(1) Suction Catheter: 8 FR
(2) IV Extension Set: Male to Male	(1) Nasogastric Tube: 5 FR - 15"
(2) 3-Way Stopcock	(1) Nasogastric Tube: 8.0 FR - 15"
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Lubrication Packet
(5) Sodium Chloride Flush Syringe: 10 mL	(1) Syringe: 35 mL
(2) Safety Needles 1.5": 22 Gauge	(1) Nasal Airway Rusch: 12 FR
(2) Safety Needles 1.5": 25 Gauge	(1) Oral Airway #0: 50 mm
(2) Luer Lock Syringe: 60 mL	(1) Oral Airway #1: 60 mm
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	



**LC017903 - Drawer 3 - Patient Weight: 9 kg - 11 kg**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 3.5 mm
(2) IV Catheter: 22 Gauge	(1) Cuffed Endotracheal Tube: 4.0 mm
(2) IV Catheter: 24 Gauge	(1) Small Stylet: 6 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 8 FR
(2) IV Extension Set: Male to Female	(1) Suction Catheter: 10 FR
(2) IV Extension Set: Male to Male	(1) Nasogastric Tube: 8.0 FR - 15"
(2) 3-Way Stopcock	(1) Lubrication Packet
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Syringe: 35 mL
(5) Sodium Chloride Flush Syringe: 10 mL	(1) LMA # 1.5
(2) Safety Needles 1.5": 22 Gauge	(1) Nasal Airway Rusch: 16 FR
(2) Safety Needles 1.5": 25 Gauge	(1) Oral Airway #1: 60 mm
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	

**LC017904 - Drawer 4 - Patient Weight: 12 kg - 14 kg**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 4.0 mm
(2) IV Catheter: 20 Gauge	(1) Cuffed Endotracheal Tube: 4.5 mm
(2) IV Catheter: 22 Gauge	(1) Small Stylet: 6 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 10 FR
(2) IV Extension Set: Male to Female	(1) Lubrication Packet
(2) IV Extension Set: Male to Male	(1) Syringe: 35 mL
(2) 3-Way Stopcock	(1) LMA # 2.0
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Oral Airway #2: 60 mm
(5) Sodium Chloride Flush Syringe: 10 mL	
(2) Safety Needles 1.5": 22 Gauge	
(2) Safety Needles 1.5": 25 Gauge	
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	



**LC017905 - Drawer 5 - Patient Weight: 15 kg - 18 kg**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 4.5 mm
(2) IV Catheter: 18 Gauge	(1) Cuffed Endotracheal Tube: 5.0 mm
(2) IV Catheter: 20 Gauge	(1) Small Stylet: 6 FR
(2) IV Catheter: 22 Gauge	(1) Suction Catheter: 10 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 12 FR
(2) IV Extension Set: Male to Female	(1) Lubrication Packet
(2) IV Extension Set: Male to Male	(1) Syringe: 35 mL
(2) 3-Way Stopcock	(1) LMA # 2.5
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Oral Airway #2: 70 mm
(5) Sodium Chloride Flush Syringe: 10 mL	
(2) Safety Needles 1.5": 22 Gauge	
(2) Safety Needles 1.5": 25 Gauge	
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	

**LC017906 - Drawer 6 - Patient Weight: 19 kg - 23 kg**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 5.0 mm
(2) IV Catheter: 18 Gauge	(1) Cuffed Endotracheal Tube: 5.5 mm
(2) IV Catheter: 20 Gauge	(1) Medium Stylet: 10 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 10 FR
(2) IV Extension Set: Male to Female	(1) Suction Catheter: 12 FR
(2) IV Extension Set: Male to Male	(1) Lubrication Packet
(2) 3-Way Stopcock	(1) Catheter Tip Syringe: 60 mL
(5) Sodium Chloride Flush Syringe: 3 mL	(1) LMA # 3.0
(5) Sodium Chloride Flush Syringe: 10 mL	(1) Oral Airway #2: 70 mm
(2) Safety Needles 1.5": 22 Gauge	(1) Oral Airway #3: 80 mm
(2) Safety Needles 1.5": 25 Gauge	
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	



**LC017907 - Drawer 7 - Patient Weight: 24 kg - 50 kg**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 5.5 mm
(2) IV Catheter: 18 Gauge	(1) Cuffed Endotracheal Tube: 6.0 mm
(2) IV Catheter: 20 Gauge	(1) Medium Stylet: 10 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 12 FR
(2) IV Extension Set: Male to Female	(1) Suction Catheter: 14 FR
(2) IV Extension Set: Male to Male	(1) Lubrication Packet
(2) 3-Way Stopcock	(1) Catheter Tip Syringe: 60 mL
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Nasal Airway Rusch: 26 FR
(5) Sodium Chloride Flush Syringe: 10 mL	(1) Nasal Airway Rusch: 28 FR
(2) Safety Needles 1.5": 22 Gauge	(1) Oral Airway #2: 70 mm
(2) Safety Needles 1.5": 25 Gauge	(1) Oral Airway #3: 80 mm
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	

**LC017908 - Drawer 8 - Patient Weight: 50 kg +**

VASCULAR ACCESSORIES	AIRWAY EQUIPMENT
(1) IV Start Kit	(1) Cuffed Endotracheal Tube: 7.0 mm
(2) IV Catheter: 14 Gauge	(1) Cuffed Endotracheal Tube: 8.0 mm
(2) IV Catheter: 18 Gauge	(1) Large Stylet: 12 FR
(2) IV Catheter: 20 Gauge	(1) Suction Catheter: 12 FR
(2) Butterfly Needle: 23 Gauge	(1) Suction Catheter: 14 FR
(2) IV Extension Set: Male to Female	(1) Lubrication Packet
(2) IV Extension Set: Male to Male	(1) Catheter Tip Syringe: 60 mL
(2) 3-Way Stopcock	(1) LMA # 4
(5) Sodium Chloride Flush Syringe: 3 mL	(1) Nasal Airway Rusch: 30 FR
(5) Sodium Chloride Flush Syringe: 10 mL	(1) Nasal Airway Rusch: 32 FR
(2) Safety Needles 1.5": 22 Gauge	(1) Oral Airway #4: 90 mm
(2) Safety Needles 1.5": 25 Gauge	(1) Oral Airway #5: 100 mm
(2) Luer Lock Syringe: 60 mL	
(2) Luer Lock Syringe: 20 mL	
(2) Luer Lock Syringe: 10 mL	
(2) Luer Lock Syringe: 3 mL	
(2) Luer Lock Syringe: 1 mL	